

RUGGED FEATURES

with Maximum Versatility





Delta Tiller

The **Pursuit**® offers a rugged design that is extremely versatile. Enjoy an adjustable delta tiller with wraparound handles for easy, one-handed operation. Equipped with front and rear suspension and large 13" drive tires, the Pursuit navigates outdoor terrain with ease. Cruise at speeds up to 9.5 mph (with HS version) and experience a full lighting package, including turn signals and an angle-adjustable LED headlight, for optimal illumination.



4-wheel Model # S713ES CA
4-wheel Model # S713ES HS

Specifications

Weight Capacity	350 lbs.
Maximum Speed ¹	Up to 6 mph HS: Up to 9.5 mph
Ground Clearance ²	3.5"
Turning Radius	68.5"
Overall Length ²	53.5"
Overall Width ²	25"
Seat-to-Deck Height	17.75"-18.75"
Seat-to-Ground Height	24.75"-26.25"
Front Tires	13" pneumatic (standard)
Rear Tires	13" pneumatic (standard)
Suspension	Front and rear
Range Per Charge ¹	Up to 33 miles HS: Up to 22 miles
Total Weight Without Batteries	161 lbs.
Heaviest Piece When Disassembled	N/A
Standard Seating	Type: Reclining high-back Weight: 43 lbs. Material: Black Vinyl Dimensions: 18" width, 18" depth
Drivetrain	Rear-wheel drive, sealed transaxle, 24-volt DC motor
Brakes	Regenerative and electromechanical
Battery Requirements ^{3,4}	Size: (2) Group 24, 75AH Weight: 53 lbs. each
Battery Charger	Off-board, 8 amp
Colors	 

Red Blue

Features

- Easy-to-read LED battery meter
- Reclining high-back seat with sliders, swivel, height adjustable headrest, flip-up, height and width adjustable armrests
- Delta tiller with wraparound handles
- Easily accessible tie-down points (for transport of unoccupied scooter)
- Easy-access tiller-mounted charger port
- Front and rear suspension
- Full lighting package includes directional signals and large lower headlight for pathway illumination
- Infinitely-adjustable tiller with conveniently located handle
- Large 13" pneumatic tires
- Powerful drivetrain; features 140 amp microprocessor-based controller
- 24-volt, 4-pole DC motor for increased power
- Jumbo front basket (standard)

Accessories

- Cell phone holder
- Cup holder
- 50" Lap belt
- Rear basket
- Rear view mirror
- Saddlebag
- Safety flag
- Walker holder
- Weather cover
- XLR USB charger

¹ Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).

² Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).

³ AGM or gel cell type required.

⁴ Battery weight may vary based on manufacturer.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Dimensions listed are for scooter only. Overall measurements will vary based on seating and accessory selections.