

Contemporary Features ENHANCE your Active Lifestyle



4-Wheel shown
in Candy Apple Red



3-Wheel shown
in **New** Ocean Blue

VICTORY® 10.2

The **Victory® 10.2** series offers great convenience and versatility with exclusive feather-touch disassembly and a delta tiller with wrap-around handles. Contemporary features such as the integrated storage/cup holder on the tiller ensure ample storage options, while a new under-tiller puddle light offers improved visibility to both the storage area and USB charging port. A new ambient voltmeter adjusts brightness to suit the user's environment, and frosted LED headlights and a low-mounted LED pathway light provide plenty of illumination in low-light environments.



PRIDE®
MOBILITY
SCOOTERS

| pridemobility.ca |





Specifications 3W: S6102 4W: S7102

Weight Capacity	400 lbs.
Maximum Speed ¹	Up to 5.2 mph
Ground Clearance ³	2.4" at transaxle, 2.1" at anti-tip, 3" at side deck
Turning Radius ³	3W: 42" 4W: 52.25"
Overall Length ³	3W: 45.25" 4W: 47"
Overall Width ³	22.25"
Seat-to-Ground Height	22" – 25"
Seat-to-Deck Height	15.25" – 18.25"
Front Tires	3.5" x 10" solid
Rear Tires	3.5" x 10" solid
Suspension	N/A
Range Per Charge ^{1,2}	3W/4W: Up to 11 miles (U-1) 3W: Up to 13 miles (40 AH) 4W: Up to 13.5 miles (40 AH)
Total Weight Without Batteries ^{3,7}	3W: 130.25 lbs. 4W: 139 lbs.
Heaviest Piece When Disassembled ³	3W: 56.5 lbs. (front section) 4W: 64.5 lbs. (front section)
Standard Seating	Type: Foldable molded plastic Weight: 28.5 lbs. Material: Black Vinyl Dimensions: 18" width, 17" depth or 20" width, 18" depth
Drive System	Rear-wheel drive, sealed transaxle with 24-volt DC motor
Dual Braking System	Regenerative and electromechanical
Battery Requirements ^{3,4,5}	Size: (2) U-1 standard; Weight: 24.5 lbs. each Size: (2) 40 AH (optional) Weight: 32.5 lbs. each
Battery Charger	3.5-amp, off-board (U-1), 5-amp, off-board (40 AH)
Colors	 

Candy Apple Red Ocean Blue

¹ Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tire condition. This specification can be subject to a variance of +10%, -5%.
² Tested in accordance with ANSI/RESNA, WC Vol 2, Section 4 & ISO 7176-4 standards. Results derived from theoretical calculation based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
³ Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.
⁴ AGM or gel cell type required.
⁵ Battery weight may vary based on battery manufacturer.
⁶ Includes standard seat weight.

Features

- A standard front basket provides ample storage
- Feather-touch disassembly permits simple frame separation with only one hand for easy transport and storage
- Charger port located on the tiller enables you to conveniently charge your scooter
- Delta tiller with ergonomic wraparound handles enables easy-to-use left or right hand operation and the ability to rest your wrist on the handle for added comfort
- USB charging port in the tiller lets you power your smart devices while on the go
- **New** under-tiller puddle light improves visibility to USB charging port and tiller-mounted storage
- **New** integrated storage/cup holder in tiller
- **New** ambient voltmeter adjusts brightness to suit your environment
- **New** frosted LED headlights and low-mounted LED pathway light brightens your route
- Pinchless hinge seat with swivel and sliders for increased comfort. Available in sizes 18" x 17" and 20" x 18"
- Easy-grip tiller adjustment knob
- Easy access tiller-mounted charger port

Accessories

- Cane clip attachment
- Cane/crutch holder (single, double)
- Cell phone holder
- Cup holder
- Forearm crutch holder
- 50", 60", 70" Lap belt
- Military patches
- Oxygen tank holder
- Quad cane holder
- Rear basket
- Rearview mirror
- Saddlebag
- Safety flag
- Walker holder
- Weather cover
- Wishbone crutch holder

Options

- 40 AH batteries with 5-amp, off board charger

Optional Seating

- Limited recline, black vinyl, high back with height adjustable headrest, flip-up height and width adjustable armrests
 • 18" x 18" (sliders) • 20" x 20" (no sliders) • 22" x 18" (no sliders)
- Deluxe contour high back black vinyl with height adjustable headrest, flip-up height and width adjustable armrests (no sliders)
 • 18" x 18" • 18" x 20" • 20" x 20" • 22" x 20"

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Dimensions listed are for scooter only. Overall measurements will vary based on seating and accessory selections.

