

## 4-Wheel

Model Number	S710LXW
Weight Capacity	400 lbs.
Maximum Speed <sup>1</sup>	Up to 8 mph
Ground Clearance <sup>3</sup>	2.2" at motor, 3.15" at mid-frame
Turning Radius <sup>3</sup>	61.25"
Overall Length <sup>3</sup>	47"
Overall Width <sup>3</sup>	25.62"
Seat-to-Ground Height	23" - 25"
Seat-to-Deck Height	16.5" -18.5"
Front Tires	3.5" x 10" solid
Rear Tires	3.5" x 10" solid
Suspension	Front and rear
Range Per Charge <sup>1,2</sup>	Up to 19.7 miles at 200 lbs. Up to 13 miles at 400 lbs.
Total Weight Without Batteries <sup>5</sup>	160 lbs.
Heaviest Piece When Disassembled <sup>3</sup>	77 lbs. (front section)
Standard Seating	Type: Reclining high-back Weight: 43 lbs. Material: Black Vinyl Dimensions: 18" width, 18" depth
Drive System	Rear-wheel drive, sealed transaxle with 24-volt DC motor
Dual Braking System	Regenerative and electromechanical
Battery Requirements <sup>4,6</sup>	Size: (2) 40 AH (Interceptor Recommended) Weight: 30.2 lbs. each
Battery Charger	5 amp, off-board
Warranty	Frame: lifetime limited Drivetrain: 2-year limited Electronics: 2-year limited Battery: 1-year limited

<sup>1</sup> Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tire condition. This specification can be subject to a variance of (+ or -) 10%.

<sup>2</sup> Tested in accordance with ANSI/RESNA, MC Vol2, Section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.

<sup>3</sup> Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.

<sup>4</sup> Battery weight may vary based on battery manufacturer.

<sup>5</sup> Includes standard seat weight.

<sup>6</sup> AGM or gel cell type required.

**PRIDE**<sup>®</sup>  
MOBILITY  
LIVE YOUR BEST<sup>®</sup>

Beamsville, Ontario • 888-570-1113

pridemobility.ca |    

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.

VICLXSPORT\_020119

**VICTORY<sup>®</sup>LX**  
**SPORT**

ENJOY A  
**SMOOTH**  
AND  
**SPORTY**  
**RIDE**



*\*Please note blue model is no longer available.*

**LIVE YOUR BEST<sup>®</sup>**

pridemobility.ca |    

# SMOOTH, COMFORTABLE RIDE

*\*Please note blue model is no longer available.*

18" x 18" Vinyl reclining high-back seat with sliders, height adjustable headrest and flip-up, height and width-adjustable armrests

Puddle light under tiller

Standard front basket

Full LED light with headlight, rear lights, brake lights and directional signals

Easy-to-grip tiller adjustment knob

Storage compartment under the tiller

CTS Suspension includes adjustable shocks for greater comfort



Automatic Brake Lights



Tiller Console



CTS Suspension

Maintain your active lifestyle with the **Victory® LX Sport**. Equipped with **Comfort-Trac Suspension (CTS)**, the Victory LX Sport easily handles more varied terrain while providing a smooth and comfortable ride. Enjoy convenient storage and a courtesy light under the tiller. A long-lasting LED headlight ensures optimal pathway illumination. Experience great features and a sporty style with the Victory LX Sport.

- CTS Suspension includes adjustable shocks for greater comfort
- Exclusive one-hand feather touch disassembly
- Charger port located on the tiller lets you conveniently charge your scooter
- Delta tiller with ergonomic wraparound handles lets you operate the scooter with one hand and rest your wrist
- Easy-to-read LED voltmeter
- Rearview mirror

## Available Colors



Candy Apple Red

## Optional Seating

Limited recline, black vinyl, high back with height adjustable headrest, flip-up height and width adjustable armrests

- 22" x 18"
- 20" x 20"

## Accessories

- Cane/crutch holder (single, double)
- Cell phone holder
- Cup holder
- Forearm crutch holder
- Oxygen tank holder
- Quad cane holder
- Rear basket
- Additional rearview mirror
- Saddlebag
- Safety flag
- Walker holder
- Weather cover
- Wishbone crutch holder
- XLR USB charger
- Lap belts: 50", 60" and 70"